

Vanderbilt Group Introduction

Julia Velkovska, Vicki Greene, Sourav Tarafdar

EIC Tracking R&D meeting

09/21/2020




Vanderbilt group composition

- Two Professors (Julia Velkovska and Vicki Greene)
- One research assistant professor (Sourav Tarafdar) and two postdocs
- 6 graduate students
- 1-2 undergraduates (more can be recruited as needed)

Research


- **PHENIX :**
 - Construction and operation of Pad-Chambers and Time of Flight West
 - Identified particle production and flow measurements, including flow of heavy-flavor particles
- **CMS :**
 - Tier 2 computing
 - Measurements of collective flow, strange particles, high-pt charged particles, jets, jet substructure
- **sPHENIX :**
 - sPHENIX TPC GEM assembly factory
 - Tracking software
- **JetScape/X-scape collaboration**
 - Model-to-data comparisons; studies of collectivity and jet-quenching in small systems

Vanderbilt MPGD Research Facility



Available equipment (not all listed) :

1. MPGD R&D bench
2. 8+4 channel two CAEN HVPS
3. 10X10 cm² multilayer GEM detector
4. APV cards with SRS
5. Mini X-ray tube for IBF measurement
6. DAQ (scope + SRS)
7. 10 channel Zagreb picoammeter
8. In house build 3 gas mixing unit (mixing inaccuracy of +/- 0.5% as on right hand plot)




Thanks to Bob Azmoun for suggestions

Additional engineering support : Availability of Physics department machine shop


Results from multilayered GEM detector in Vanderbilt MPGD lab

- Several studies for triple and quad GEM detectors were done for ArCO₂ gas mixtures using 10x10 standard GEMs.
- Initial study used voltage divider for biasing detector electrodes.
- Spectra is measured from bottom of last GEM using electronics chain of preamp + shaping amp + scope

Results 1 : Using voltage divider




MDO3054 - 4:09:47 PM 2/10/2020






MDO3054 - 10:20:40 PM 6/13/2020

Results 2 : Using individual channel of HVPS on Quad GEM




MDO3054 - 4:57:22 PM 9/18/2020



IBF in Quad gem for ArCO₂(70:30) using X-ray tube.

- Anode and cathode current
- IBF and effective gain correlation



Activities on track :

- In the process of getting Fe-55 spectra using SRS (Thanks to Martin for setting it up !).
- Study the effect of gap fields on IBF .
- Using misaligned gem holes for IBF blocking.
- Performance study of hybrid MPGD

Future Plans

- **The Vanderbilt MPGD lab is well equipped for performing MPGD related R&D for EIC tracking**
- **Hardware and software expertise, and track record of completing challenging projects**
- **The group is committed to EIC research, including detector R&D, construction, and eventually – data analysis**
- **Close collaboration with BNL (Craig Woody et.al.) is envisioned for studying different readout pads for TPC using GEMs, and also a hybrid MPGD detector.**